

Paediatric Normal Values - Anaesthesia

All codes listed are Intersurgical product codes

Interact with us

www.intersurgical.co.uk



Age	Weight (kg)	Height or length (cm)	Anaesthetic facemask size	Oropharyngeal airway size (ISO) *	i-gel size	LM size	LM cuff maximum inflation value (ml)	Tracheal tube uncuffed internal diameter (mm) **	Tracheal tube cuffed internal diameter (mm)	Awake heart rate (per minute)	Systolic blood pressure (mm Hg)	Respiratory rate (per minute)	Estimated tidal volume (ml)	Approximate blood volume (ml)	Acceptable haematocrit	Intravenous maintenance fluid (ml/hr) ***	Breathing system HMEF	Anaesthetic circuit		
< 1 month	3.5	50	0 1127000	00 (5.0) 1110050	1 8201000	1 8001000	< 4	3.0		85 - 205	60	30 - 40	21	315	≥ 0.30	14				
1 month	4	54						3.0 - 3.5	3.0		80		24						360	16
3 months	5	60	1 1128000	0 (5.5) 1110055	1.5 8215000	1.5 8015000	< 7	3.5	3.0	100 - 180	80	30 - 40	30	400	≥ 0.25	20	Clear-Therm Micro (1441000)	Jackson-Rees T-piece with 0.5 litre reservoir bag (2121000, 2122000 - with APL valve) in anaesthetic room and 15 mm paediatric circle system (2142000) in theatre		
6 months	7	67						4.0	3.5		80		42						560	28
1 year	10	78	2 1129000	1 (6.5) 1111065	2 8202000	2 8002000	< 10	4.0	3.5	100 - 180	92	24 - 30	60	800	≥ 0.2	40	Clear-Therm Mini (1831000)			
2 years	12	87						4.5	4.0		94		72						900	44
3 years	14	95						4.5 - 5.0	4.0 - 4.5		96		84						1050	48
4 years	16	103						5.0	4.5		98		96						1200	52
5 years	18	109						5.0 - 5.5	4.5 - 5.0		100		108						1350	56
6 years	20	116	3 1123000	1.5 (7.0) 1111570	2 8202000	2.5 8025000	< 14	5.5	5.0	60 - 140	102	20 - 24	120	1500	≥ 0.2	60	Clear-Therm Mini (1831000)	Mapleson C in anaesthetic room (2108000 or 2102000), adult circle system in theatre (2010000)		
7 years	22	122						5.5 - 6.0	5.0 - 5.5		104		132						1650	62
8 years	26	128						6.0 - 6.5			106		156						1820	66
10 years	30	139	4 1124000	2 (8.0) 1112080	2.5 8225000	3 8003000	< 20	7.0		60 - 100	110	20 - 24	180	2100	≥ 0.2	70	Clear-Therm 3 (1541000)			
12 years	38	149						7.0 - 7.5			114		228						2660	78
Adolescent	50	161	3 8203000	4 8004000	< 30	7 - 8	118		12 - 20	118	300	3500	90							
Reference	1	2					3	3		3	3	3	1	1	1	1	1	4	4	6

Ambient temperature to be a minimum of 21 degrees Celsius. For smaller children and neonates undergoing surgery or resuscitation, additional warming with a Bair Hugger® or similar device is recommended (5).
 * The correct size of an oropharyngeal airway is one that, when laid against the side of the face, has a length equal to the distance between the level of the patient's incisors (or where they will be) to the angle of the jaw (1).
 ** A correctly sized uncuffed tracheal tube should have a small audible leak around the tube when 20 cm of water pressure is applied from the breathing system (5).

*** Intravenous maintenance fluid recommendations for previously well children aged from one month to 16 years old

The majority of children may be safely administered sodium chloride 0.45% with glucose (2.5 or 5%). Do not use sodium chloride 0.18% with glucose 4%.
 Some children at high risk of hyponatraemia should only receive isotonic fluids (see list opposite).
 Some acutely ill children with increased anti-diuretic hormone (ADH) secretion (e.g. post-operative patients or those with intracranial infections or head injuries) may benefit from their maintenance fluid being restricted to two-thirds normal recommended volume.
 To avoid dangerous hypo or hypernatraemia, monitor the child's weight and calculate fluid balance. Use a volumetric pump. Check plasma electrolyte and glucose concentration before and regularly throughout intravenous therapy.
 Consider adding potassium 40 mmol/l to maintenance fluids once plasma potassium levels are known.
 Children requiring both maintenance fluids and replacement of ongoing losses should receive a single isotonic fluid.

Children who should only receive isotonic fluids include those who:

are peri- or post-operative
 have low plasma sodium
 have CNS infection or a head injury
 have sepsis
 have a self-wasting syndromes
 require the replacement of ongoing losses
 have intravascular volume depletion or hypotension
 have bronchiolitis
 have excessive gastrointestinal losses
 have a chronic condition such as diabetes, cystic fibrosis or a pituitary deficit
Examples of isotonic fluids are: sodium chloride 0.9%, sodium chloride 0.9% with 5% glucose or Hartmann's solution.
For further information regarding the treatment of shock and the replacement of pre-existing fluid deficit, consult the NPSA website, EPALS manual and other appropriate resources.

References

1. European Paediatric Life Support, 4th Edition, 2016; p3-4, 37, 149. Reproduced with the kind permission of the Resuscitation Council (UK).
2. Child Growth Foundation (Charity Reg No 274325). Boys and Girls Growth Charts (Birth-18 years). London: 2 Mayfield Avenue, London. W4 1PW, 1996.
3. Intersurgical product information <http://www.intersurgical.co.uk/>. Accessed June 2017.
4. Cunliffe M. Fluid and electrolyte management in children. BJA CEPD reviews 2003; 3(1): p1-4.
5. Basic techniques for anaesthesia. In Sumner E and Hatch DJ, eds. Paediatric Anaesthesia. London: Arnold, a member of the Hodder Headline Group 2000; p182, 194.
6. Reducing the risk of hyponatraemia when administering intravenous infusions to children, March 2007; p11. <http://www.nrls.npsa.nhs.uk/resources/?EntryId45=59809>. Accessed June 2017.

Disclaimer

Whilst every care has been taken to ensure that doses and recommendations are correct, the responsibility for final checking must rest with the practitioner. The authors cannot accept any responsibility for errors in this publication. Equipment sizes are based upon Intersurgical recommendations.
 Copyright
 © Intersurgical 2018, all rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior written permission of the copyright holder.
 Accreditation
 Intersurgical would like to thank Richard Hixson, Consultant Anaesthetist at Darlington Memorial Hospital for the concept and his help in the production of this chart.