



PEEP valves and accessories for the StarMed hoods



Critical Care ▪ Patient Interfaces ▪ **STARMed**

PEEP valves

Adjustable PEEP valves for CPAP therapy



| Code | Description | Box Qty. |
|---------|-----------------------|----------|
| SPV20 | Yellow and blue color | 10 |
| 2226000 | Green and blue color | 10 |

Accessories



| Code | Description | Box Qty. |
|---------|---|----------|
| STR/028 | Nebuliser accessory, complete with mouthpiece for delivering nebulised drugs during CPAP or NIV therapy | 5 |
| STR/014 | 22F/22F straight adaptor for connecting the hood to some breathing systems configurations | 10 |
| STR/019 | 22F cap to occlude the expiratory side of the hood in some breathing systems configurations | 10 |
| CP903 | Cuff inflation bulb for inflating the built in neck cushion present in some hood ranges | 5 |
| CP908 | Earplugs for improving patient comfort during therapy | 25 |
| CP917 | 22F/22F adaptor with built in manometer to convert any NIV hood into a CPAP hood | 5 |
| SPV108 | 22F/22F adaptor with built in manometer and adjustable PEEP valve to convert any NIV hood for immediate use in CPAP therapy | 10 |
| CP900 | Spare underarm straps. Spare standard fixations for adult CaStar and Ventukit hoods | 5 |
| CP920 | Spare underarm straps for CaStar Up. Spare fixations for adult CaStar Up openable hood series | 5 |
| CP906 | Elastic belts with slots for bed fixation. An alternative hood fixation system | 5 |
| CP912 | Infant harness – Low. Spare fixation for the Infant hood, smaller patient size | 5 |
| CP913 | Infant harness – High. Spare fixation for the Infant hood, larger patient size | 5 |



www.intersurgical.co.uk/info/starmed

Make an enquiry



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



Manufacturer: Intersurgical S.p.A. - Mirandola (MO) - Italy
 The Intersurgical S.p.A. Organisation has a Quality Management System Certificate by
 KIWA CERMET ITALY S.p.A. according to ISO 9001:2015 and ISO 13485:2016

Please think before you print
 Save energy and paper. If you must print this
 information sheet please print it double sided.

IS1.47_INT • Issue 2 02.19