



# ClearLite™ anaesthetic face masks

A clear vision for the future 




**Anaesthesia** ▪ Anaesthetic Face Masks

## ClearLite™ anaesthetic face masks

A range of anaesthetic face masks which includes seven sizes to meet the clinical needs of all patient types and sizes. The range is available in sizes 0 neonatal to size 6 large adult. Due to the enhanced anatomical design, the ClearLite mask offers superior seal performance with minimum pressure, reducing the risk of unnecessary discomfort or potential trauma. The superior seal will also assist the clinician in

providing optimum manual ventilation of the patient in both the anaesthetic and emergency situation.

ClearLite is non-PVC and phthalate free which significantly minimises the impact upon the environment. 

### Lightweight mask

Provides enhanced grip and comfort in use for the clinician

### Anatomical shape

Provides a precise fit for different face types and sizes. Minimum pressure required to achieve an effective seal reducing the risk of inadvertent trauma and nerve damage

### Crystal clear body

Provides excellent visibility for confirmation that the patient is well perfused

### Reduced mask volume

Reduces potential rebreathing and risks associated with elevated levels of end tidal CO<sub>2</sub>



### Safe and secure connection

This is achieved between the mask and the breathing system in use. The ISO 22mm tapered connection complies with ISO 5356

### Reduce waste

Available without a hook ring

### Flexible non-inflatable cuff

Provides a comfortable and effective airtight seal


### Colour-coding

The comfortable, flexible seal is colour coded to allow quick and easy identification of the correct sized mask















### Low environmental impact

Non-PVC and phthalate free



 Video available at [www.intersurgical.com](http://www.intersurgical.com)



Code	Code without hook ring	Description	Box Qty.
 7296000	7296001	ClearLite, anaesthetic face mask, size 6, extra large adult, red seal, 22F	20 
 7295000	7295001	ClearLite, anaesthetic face mask, size 5, large adult, orange seal, 22F	35 
 7294000	7294001	ClearLite, anaesthetic face mask, size 4, adult, green seal, 22F	35 
 7293000	7293001	ClearLite, anaesthetic face mask, size 3, small adult, yellow seal, 22F	35 
 7292000	7292001	ClearLite, anaesthetic face mask, size 2, paediatric, white seal, 22F	25 
 -	7291000	ClearLite, anaesthetic face mask, size 1, infant, grey seal, no hook ring, 15M	30 
 -	7290000	ClearLite, anaesthetic face mask, size 0, neonate, blue seal, no hook ring, 15M	40 

Make an enquiry



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 [info@intersurgical.com](mailto:info@intersurgical.com) [www.intersurgical.com](http://www.intersurgical.com)



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.

IS5.31 • Issue 2 07.19