Intersurgical Guedel Airways

New ISO 5364:2016 Specified Colours

| Product Code | Intersurgical Size | ISO Size | Previous Intersurgical Colour | | New ISO Compliant Colour | |
|--------------|-----------------------|----------|----------------------------------|--|-----------------------------|-----|
| 1000035 (S*) | 000 | 3.5 | Pink | | Light Green | NEW |
| 1100050 (S*) | 00 | 5.0 | Blue | | Blue | |
| 1110055 (S*) | 0 | 5.5 | Grey | | Grey | |
| 1111065 (S*) | 1 | 6.5 | White | | Brown | NEW |
| 1111570 (S*) | 1.5 | 7.0 | Yellow | | White | NEW |
| 1112080 (S*) | 2 | 8.0 | Green | | Green | |
| 1113090 (S*) | 3 | 9.0 | Orange | | Yellow | NEW |
| 1114100 (S*) | 4 | 10.0 | Red | | Red | |
| 1115120 (S*) | 5 | 12.0 | Purple | | Purple | |

This colour change effects all of Intersurgical's Guedel Airway products, including bulk-packed options.

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S

Colour change requirement

In order to comply with the colour coding system specified in ISO 5364:2016, it was necessary for Intersurgical to update four sizes of our Guedels. The sizes that have changed are as follows:

- ISO size 3.5; from Pink to Light Green
- ISO size 6.5; from White to Brown
- ISO size 7.0; from Yellow to White
- ISO size 9.0; from Orange to Yellow

Each department will need to consider locally how this change will be implemented, as it may take time for existing stock to be used.









The four sizes of Intersurgical Guedel airways that have had their colours changed to conform to thea updated standard.

For more information, please visit www.intersurgical.com/info/guedelairway



