



Right-angled water traps

Self-sealing water traps



Critical Care ▪ Accessories



Quality, innovation and choice

Interact with us



www.intersurgical.com

Right-angled water traps

Self-sealing water traps

Intersurgical's right-angled water traps for the collection of condensation compliments our existing range of products, the angled water trap is available in two sizes.

Medium volume and large volume, both are ideal for use in applications where a 90° attachments may be required.

The self-sealing facility allow the water traps to be emptied without interrupting gas flow.

Both water traps feature a 22F machine side connection to a 22M breathing system side connection, and can be attached directly to a ventilator or to an anaesthetic machine.



Self-sealing facility

Allow the water traps to be emptied without interrupting flow



Code	Description	Box Qty.
1916000	Water trap, self-sealing, right-angled, medium volume	55
1915000	Water trap, self-sealing, right-angled, large volume	40

Make an enquiry



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
 Save energy and paper.
 If you must print this information sheet please print it double sided.