



## Coach<sup>®</sup> 2 Volume Incentive Spirometers

Makes incentive spirometry easier for both patients and staff



Oxygen and Aerosol Therapy ▪ Respiratory Physiotherapy

## Coach® 2 Volume Incentive Spirometers

The Coach® 2 volumetric incentive spirometers are designed to help patients who have undergone thoracic or abdominal surgery. The patient is encouraged to take long, slow, deep breaths. The Coach® 2 measures these breaths and visually indicates the volume achieved during sustained maximal inspiration. Incentive spirometry is intended to increase the inspiratory volume and to prevent post-operative pulmonary complications.

Models are available for adults and children.

### Volume scale printed on both sides

enables clinician to view patient's performance

### Single patient use

reduces cross contamination



### Integral one-way valve

ensures patients inhale, rather than exhale into the unit

### Instruction book in the base

help patient's with their exercises without direct supervision

### Integral connector

provides supplemental oxygen

### Adjustable volume indicator

enables targets to be easily set

Code	Description	Volume	Box Qty.
1535000	Coach® 2 Adult with one-way valve	4000ml	12
1536000	Coach® 2 for Kids with one-way valve	2000ml	12

Make an enquiry

Coach is a registered trademark of Smith Medical.

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS10.25 • Issue 3 04.18



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
 T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com

UK • Ireland • France • Germany • Spain • Portugal • Italy • Benelux • Sweden • Denmark • Lithuania • Russia • Czech Republic  
 Turkey • South Africa • China • Japan • Taiwan • Philippines • USA • Canada • Colombia • Australia